HIPPEX CRYO/HIPPEX FREEZE

Freezing Media for Equine Sperm

Operating instructions

The general guidelines about sperm freezing like for example centrifugation time and speed, target sperm concentration, straws filling, freezing methods and freezing protocols are considered as being known.

1. Place the **HIPPEX GLYCEROL Solution** bottle (50 ml) in the water bath (or incubator) at 37°C.
2. Add the **HIPPEX GLYCEROL Solution** (50 ml at 37°C) into the bottle **HIPPEX CRYO** or **HIPPEX FREEZE**.
3. After gentle shaking and when the powder is completely dissolved, keep the **HIPPEX CRYO** or **HIPPEX FREEZE** at 37°C for another 30-60 minutes to stabilize the solution. From time to time make a gentle shaking or place the solution on a roller bank. After 30 to 60 minutes, the solution is ready to be used.
4. Resuspend the pellet of the centrifuged sperm (centrifugation was done by using HIPPEX Centrifugation Medium, HCM-D1) to the necessary desired target sperm motile concentration by adding the correct volume of **HIPPEX CRYO** or **HIPPEX FREEZE**. Place the obtained mixture on a shaking table or on a roller bank for 15 minutes to optimize the mixing and to obtain equilibrium.
5. Follow afterwards the protocols that are used in your lab for filling the straws and freezing procedure.

*PS: The delivered **HIPPEX GLYCEROL Solution** contains **2.4%** glycerol. If needed, the glycerol concentration can be adapted on demand.*

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HIPPEX CRYO & FREEZE OPERATING INSTRUCTIONS – Version 4.10– Jan 2023 –

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Copyright © 2023 HUVESearch BV, Belgium